

# Cygnus EndonucleaseGTP® Assay

## Immunoenzymetric Assay for the Measurement of Endonuclease Catalog # F960

### Intended Use

This kit is intended for use in determining the presence of any endonuclease impurities in products manufactured by virus production. The kit is for **Research and Manufacturing Use Only** and is not intended for diagnostic use in humans or animals.

### Summary and Explanation

Expression of recombinant viral vectors and vaccines in mammalian cell culture is a cost-effective method for production of commercial quantities of these novel biological drugs. The manufacturing and purification process of these products leaves the potential for impurities by HCPs and other additives such as enzymes used to remove endogenous DNA and RNA, and plasmid DNA used for viral vector production. Such impurities can reduce the efficacy of the therapeutic agent and result in adverse toxic or immunological reactions, and thus it is desirable to reduce HCP and added enzyme impurities to the lowest levels practical.

This simple to use, objective, and semi-quantitative ELISA is a sensitive and specific method to aid in optimal purification process development, process control, and in routine quality control and product release testing. The antibodies have been generated against and affinity purified against endonuclease. As such, this kit can be used as a process development tool to monitor the optimal removal of endonuclease impurities as well as in routine final product release.

This highly sensitive ELISA kit developed with these antibodies is an accurate, sensitive, and precise method for detection of genetically engineered endonuclease from *Serratia marcescens* such as Benzonase® Nuclease and DENARASE® in samples throughout the purification process. Each user of this kit is encouraged to perform a similar qualification study to demonstrate it meets their analytical needs

### Principle of the Procedure

The endonuclease assay is a two-site immunoenzymetric assay. Samples containing endonuclease are reacted simultaneously with a horseradish peroxidase (HRP) enzyme labeled anti-endonuclease antibody (goat polyclonal) in microtiter strips coated with an affinity purified capture anti-endonuclease antibody (goat

polyclonal). The immunological reactions result in the formation of a sandwich complex of solid phase antibody-endonuclease-enzyme labeled antibody. The microtiter strips are washed to remove any unbound reactants. The substrate, tetramethylbenzidine (TMB) is then reacted. The amount of hydrolyzed substrate is read on a microtiter plate reader and is directly proportional to the concentration of endonuclease enzyme present.

### Reagents & Materials Provided

Component	Product #
<b>Anti-EndonucleaseGTP®:HRP</b> Affinity purified goat antibody conjugated to HRP in a protein matrix with preservative. 1x12mL	<b>F961</b>
<b>Anti-EndonucleaseGTP® coated microtiter strips</b> 12x8 well strips in a bag with desiccant	<b>F962</b>
<b>EndonucleaseGTP® Standards</b> Endonuclease in a protein matrix with preservative. Standards at 0, 0.31, 0.63, 1.25, 2.5, 5, 10, and 20ng/mL. 1 mL/vial	<b>F963</b>
<b>Stop Solution</b> 0.5M sulfuric acid. 1x12mL	<b>F006</b>
<b>TMB Substrate</b> 3,3',5,5' Tetramethylbenzidine. 1x12mL	<b>F005</b>
<b>Wash Concentrate (20X)</b> Tris buffered saline with preservative. 1x50mL	<b>F004</b>

### Storage & Stability

- All reagents should be stored at 2°C to 8°C for stability until the expiration date printed on the kit.
- Reconstituted wash solution is stable until the expiration date of the kit.
- After prolonged storage, you may notice a salt precipitate and/or yellowing of the wash concentrate. These changes will not impact assay performance. To dissolve the precipitate, mix the wash concentrate thoroughly and dilute as directed in the 'Preparation of Reagents' section. Reconstituted wash solution is stable until the expiration date of the kit.

## Materials & Equipment Required But Not Provided

- Microtiter plate reader spectrophotometer with dual wavelength capability at 450 & 650nm. *(If your plate reader does not provide dual wavelength analysis you may read at just the 450nm wavelength.)*
- Pipettors - 50µL and 100µL
- Repeating or multichannel pipettor - 100µL
- Microtiter plate rotator (400-600 rpm)
- Sample Diluent (recommended Cat # 1600)
- Distilled water
- 1 liter wash bottle for diluted wash solution

## Precautions

- For Research or Manufacturing use only.
- Stop reagent is 0.5M H<sub>2</sub>SO<sub>4</sub>. Avoid contact with eyes, skin, and clothing.
- This kit should only be used by qualified technicians.

## Preparation of Reagents

- Bring all reagents to room temperature.
- Dilute wash concentrate to 1 liter in distilled water, label with kit lot and expiration date, and store at 4°C.

## Procedural Notes

1. Complete washing of the plates to remove excess unreacted reagents is essential to good assay reproducibility and sensitivity. The manual wash procedure described below generally provides lower backgrounds, higher specific absorbance, and better precision. If duplicate CVs are poor, or if the absorbance of the '0' standard is greater than 0.300, evaluate plate washing procedure for proper performance.

2. High Dose Hook Effect or poor dilutional linearity may be observed in samples with very high concentrations of enzyme. High Dose Hook Effect is due to insufficient excess of antibody for very high concentrations of the analyte present in samples upstream in the purification process. Samples greater than 780ng/mL may give absorbances less than the 20ng/mL standard. In such cases the absorbance of the sample at all dilutions may be lower than the highest standard in the kit, however these samples will fail to show acceptable dilution linearity as evidenced by an apparent increase in dilution corrected enzyme concentration with increasing dilution. High Dose Hook and poor dilutional linearity are most likely to be encountered from samples early in the

purification process. If a hook effect is possible, samples should also be assayed diluted. If the endonuclease concentration of the undiluted sample is less than the diluted sample this may be indicative of the hook effect. Such samples should be diluted at least to the minimum required dilutions (MRDs) as established by your qualification studies using your actual final and in-process drug samples. The MRD is the first dilution at which all subsequent dilutions yield the same enzyme value within the statistical limits of assay precision. The enzyme value to be reported for such samples is the dilution corrected value at or greater than the established MRD. The diluent used should be compatible with accurate recovery. The preferred diluent is our Cat# 1600 available in 25 or 100mL bottles. This is the same material used to prepare the kit standards. As the sample is diluted in 1600, its matrix begins to approach that of the standards, thus reducing any inaccuracies caused by dilutional artifacts. Other prospective diluents must be tested for non-specific binding and recovery by using them to dilute the 20ng/mL standard, as described in the "Limitations" section below.

## Limitations

- Before relying exclusively on this assay to detect endonuclease, each laboratory should qualify that the kit antibodies and assay procedure yield acceptable specificity, accuracy, and precision. A suggested protocol for this qualification can be obtained from our Technical Services Department or our web site.
- The standards used in this assay are comprised of known concentrations of endonuclease. The antibodies used in the development of this kit demonstrated a recovery of 1:1 ratio of known concentrations of endonuclease.
- Certain sample matrices may interfere in this assay. The standards used in this kit attempt to simulate typical sample protein and matrices. However the potential exists that the product itself or other components in the sample matrix may result in either positive or negative interference in this assay. High or low pH, detergents, urea, high salt concentrations, and organic solvents are some of the known interference factors. It is advised to test all sample matrices for interference by diluting the 20ng/mL standard, 1 part to 4 parts of the matrix containing no or very low endonuclease. This diluted standard when assayed as an unknown, should give an added value in the range of 3 to 5ng/mL. Consult *Cygnus Technologies* Technical Service Department for advice on how to quantitate the assay in problematic matrices.
- Avoid the assay of samples containing sodium azide (NaN<sub>3</sub>) which will destroy the HRP activity of the conjugate and could result in the underestimation of Endonuclease enzyme levels.

## Assay Protocol

- The protocol specifies use of an approved orbital microtiter plate shaker for the immunological steps. These can be purchased from most laboratory supply companies. If you do not have such a device, it is possible to incubate the plate without shaking however it will be necessary to extend the immunological incubation step in the plate by about one hour in order to achieve comparable results to the shaking protocol. **Do not shake during the 30 minute substrate incubation step, as this may result in higher backgrounds and worse precision.**
- Bring all reagents to room temperature.
- Set-up plate spectrophotometer to read dual wavelength at 450nm for the test wavelength and ~650nm for the reference.
- Thorough washing is essential to proper performance of this assay. The manual method described in the assay protocol is preferred for best precision, sensitivity and accuracy. A more detailed discussion of this procedure can be obtained from our Technical Services Department or on our web site. In addition, a video demonstration of proper plate washing technique is available in the 'Technical Help' section of our web site.
- All standards, controls, and samples should be assayed at least in duplicate.
- Maintain a repetitive timing sequence from well to well for all assay steps to ensure that all incubation times are the same for each well.
- Make a work list for each assay to identify the location of each standard, control, and sample.
- It is recommended that your laboratory assay appropriate quality control samples in each run to ensure that all reagents and procedures are correct. **You are strongly urged to make controls in your typical sample matrix using endonuclease derived from your process. These controls can be aliquoted into single use vials and stored frozen for long-term stability.**
- If the substrate has a distinct blue color prior to assay it may have been contaminated. If the absorbance of 100µL of substrate plus 100µL of stop against a water blank is greater than 0.1 it may be necessary to obtain new substrate or the sensitivity of the assay may be compromised.
- Strips should be read within 30 minutes after adding stop solution since color will fade over time.

## Assay Protocol

- Pipette 100µL of anti-EndonucleaseGTP®:HRP (#F961) into each well.
- Pipette 50µL of standards, controls and samples into wells indicated on work list.
- Cover & incubate on orbital shaker at 400 - 600rpm for 1 hour at room temperature, 24°C ± 4°C.
- Dump contents of wells into waste. Blot and gently but firmly tap over absorbent paper to remove most of the residual liquid. Overly aggressive banging of the plate or use of vacuum aspiration devices in an attempt to remove all residual liquid is not necessary and may cause variable dissociation of antibody bound material resulting in lower ODs and worse precision. Fill wells generously to overflowing with diluted wash solution using a squirt bottle or by pipetting in ~350µL. Dump and tap again. Repeat for a total of 4 washes. Wipe off any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step. Do not allow wash solution to remain in wells for longer than a few seconds. Do not allow wells to dry before adding substrate.
- Pipette 100µL of TMB substrate (#F005).
- Incubate at room temperature for 30 minutes. **DO NOT SHAKE.**
- Pipette 100µL of Stop Solution (#F006).
- Read absorbance at 450/650nm.

## Example Data

Well #	Contents	Abs. at 450-650nm	Mean Abs.
A1	Zero Std	0.070	0.068
A2	Zero Std	0.065	
B1	0.31ng/mL	0.109	0.109
B2	0.31ng/mL	0.109	
C1	0.63ng/mL	0.165	0.164
C2	0.63ng/mL	0.162	
D1	1.25ng/mL	0.254	0.255
D2	1.25ng/mL	0.256	
E1	2.5ng/mL	0.436	0.438
E2	2.5ng/mL	0.441	
F1	5ng/mL	0.808	0.791
F2	5ng/mL	0.774	
G1	10ng/mL	1.383	1.362
G2	10ng/mL	1.341	
H1	20ng/mL	2.480	2.436
H2	20ng/mL	2.392	

## Calculation of Results

The standards may be used to construct a standard curve with values reported in ng/mL endonuclease. This data reduction may be performed through computer methods using curve-fitting routines such as point-to-point, cubic spline, or 4 parameter logistic fit. **Do not use linear regression analysis to interpolate values for samples as this may lead to significant inaccuracies!** Data may also be manually reduced by plotting the absorbance values of the standard on the y-axis versus concentration on the x-axis and drawing a smooth point-to-point line. Absorbances of samples are then interpolated from this standard curve.

## Quality Control

- Precision on duplicate samples should yield average % coefficients of variation of less than 10% for samples in the range of 0.31-20ng/mL. CVs for samples less than 0.31ng/mL may be greater than 10%.
- It is recommended that each laboratory assay appropriate quality control samples in each run to ensure that all reagents and procedures are correct.

## Performance Characteristics

*Cygnus Technologies* has qualified this assay by Harmonization analytical criteria. Any new sample types must be qualified by your lab to determine MRD and acceptable spike & recovery as described above and in our Qualification Summary report. Operators should refer to that report for specifics on methods used in qualification, and expected assay performance. This qualification is generic in nature and is intended to supplement but not replace a comprehensive user and sample type qualification that should be performed by each laboratory.

## Sensitivity

The lower limit of detection (**LOD**) is defined as that concentration corresponding to a signal two standard deviations above the mean of the zero standard. LOD is ~0.06 ng/mL.

The lower limit of quantitation (**LLOQ**) is defined as the lowest concentration, where concentration coefficients of variation (CVs) are less than 20%. The LLOQ is ~0.31ng/mL.

## Specificity/Cross-Reactivity

Cross reactivity to non-endonuclease components has not been extensively investigated with this kit. You should evaluate components in your samples for positive interferences such as cross reactivity and non-specific binding. Negative interference studies are described below.

## Precision

Both intra (n=20 replicates) and inter-assay (n=10 assays) precision were determined on 3 pools with low (~0.75ng/mL), medium (~3ng/mL), and high concentrations (~12ng/mL). The % CV is the standard deviation divided by the mean and multiplied by 100.

Pool	Intra assay CV	Inter assay CV
Low	5.1%	2.9%
Medium	3.6%	2.1%
High	3.6%	4.3%

## Recovery/Interference Studies

Controls with known concentrations of endonuclease enzyme were used to evaluate the performance of this assay. All of these samples yielded acceptable recovery defined as between 80-120%. The standards used in this kit contain 4mg/mL of bovine serum albumin intended to simulate non-specific protein effects of most sample proteins. However, very high concentrations of some products may interfere in the accurate measurement of endonuclease analyte. In general, extremes in pH (less than 5.0 and greater than 8.5), high salt concentration, high polysaccharide concentrations, urea, organic solvents, and most detergents can cause under-recovery. Each user should qualify that their sample matrices yield accurate recovery. Such an experiment can be performed, by diluting the 20ng/mL standard provided with this kit, into the sample matrix in question as described in the "Limitations" section. *Cygnus Technologies* offers a more concentrated form of the HCP used to prepare the kits standards for your spike recovery and preparation of analyte controls.

Due to the differences in how manufacturers of Endonuclease enzymes report their enzyme specific activities in units/mg and the fact that specific activity reflects active enzymatic activity per mg protein while not accounting for the inactive enzyme mass, you can expect up to 3X recovery value for a given standard provided in the EndonucleaseGTP ELISA Kit, Item Number F960, as compared to Merck KGaA Benzonase ELISA Kit II. *Cygnus* recommends preparing a positive control using your specific endonuclease to better assess the expected differences.

## Hook Capacity

Increasing concentrations of endonuclease enzyme greater than 20ng/mL were assayed as unknowns. The hook capacity, defined as that concentration yielding an absorbance reading less than the 20ng/mL standard was ~780ng/mL.

## Ordering Information/ Customer Service

*Cygnus Technologies* also offers kits for the detection of HCPs from HEK 293 cells, as well as extraction and detection of Human Host Cell DNA. The following kits are available:

- Cat# F650S, HEK 293 HCP ELISA
- Residual Host Cell DNA extraction:  
Cat # D100W, DNA Extraction Kit in 96 deep well plate  
Cat # D100T, DNA Extraction Kit in microfuge tubes
- Residual Human Host Cell DNA extraction and detection using PicoGreen® dye:  
Cat # D160W, DNA Extraction Kit in 96 deep well plate  
Cat # D160T, DNA Extraction Kit in microfuge tubes

To place an order or to obtain additional product information contact *Cygnus Technologies*:

[www.cygnustechnologies.com](http://www.cygnustechnologies.com)  
Cygnus Technologies, LLC  
1523 Olde Waterford Way  
Leland, NC 28451 USA  
Tel: 910-454-9442

Email for all Order inquiries:  
[orders@cygnustechnologies.com](mailto:orders@cygnustechnologies.com)

Email for Technical Support:  
[techsupport@cygnustechnologies.com](mailto:techsupport@cygnustechnologies.com)

---

Benzonase® Nuclease is a registered trademark of Merck KGaA  
DENARASE® is a registered trademark of c-LEcta

