

# SF9 Insect Cell Host Cell Proteins 2<sup>nd</sup> Generation

## Immunoenzymetric Assay for the Measurement of SF9 Insect Cell Host Cell Proteins Catalog # F840

### Intended Use

This kit is intended for use in determining the presence of host cell protein contamination in products manufactured by expression in SF9 insect host cells. The kit is for **Research and Manufacturing Use Only** and is not intended for diagnostic use in humans or animals. Users should validate this assay for use with their product samples. This kit is a 2<sup>nd</sup> Generation method replacing our earlier Catalog # F020, SF9 Insect Cell HCP ELISA kit.

### Summary and Explanation

Recombinant expression by SF9 insect cells is a widely used procedure to obtain sufficient and cost effective quantities of a desired protein or virus. Many of these products are intended for use as therapeutic agents in humans and animals and as such must be highly purified. The manufacturing and purification process of these products leaves the potential for contamination by host cell proteins from SF9 insect cells. Such contamination can result in adverse toxic or immunological reactions and thus it is desirable to reduce host cell contamination to the lowest levels practical. Immunological methods using antibodies to HCPs such as Western Blot and ELISA are conventionally accepted. While Western Blot may be able to detect HCPs in samples from upstream in the purification process, it often lacks adequate sensitivity and specificity to detect HCPs in purified downstream and final product. The microtiter plate immunoenzymetric assay (ELISA) method employed in this kit overcomes the limitations of Western blots providing on the order of 100 fold better sensitivity. This simple to use, highly sensitive, objective, and semi-quantitative ELISA is a powerful method to aid in optimal purification process development, process control, routine quality control, and product release testing. This kit is "generic" in the sense that it is intended to react with essentially all of the HCPs that could contaminate the product independent of the purification process. The antibodies have been generated against and affinity purified using a mild lysate of SF9 cells to obtain HCPs typically encountered in your initial product recovery step. Coverage of the antibody to individual HCPs was determined by a method termed

Antibody Affinity Extraction (AAE). AAE is much more sensitive and specific than 2D Western blot. To determine coverage to your process specific HCPs, it is recommended to use AAE. For more information on AAE analysis please contact our Technical Services department.

Special procedures were utilized in the generation of these antibodies to insure that low molecular weight and less immunogenic contaminants as well as high molecular weight components would be represented. As such this kit can be used as a process development tool to monitor the optimal removal of host cell contaminants as well as in routine final product release testing. When the kit can be satisfactorily validated for your samples, the application of a more process specific assay is probably not necessary in that such an assay would only provide information redundant to this generic assay. However, if your validation studies indicate the antibodies in this kit are not sufficiently reactive with your process specific HCPs it may be desirable to also develop a more process specific ELISA. A process specific assay may require the use of a more specific and defined antisera. Alternatively, if the polyclonal antibody used in this kit provides sufficient sensitivity and broad antigen reactivity, it may be possible to substitute the standards used in this kit for ones made from the contaminants that typically co-purify through your purification process and thus achieve better accuracy for process specific HCPs. The suitability of this kit for a given sample type and product must be determined and validated experimentally by each laboratory. The use of a process specific assay with more defined antigens and antibodies in theory may yield better sensitivity however, such an assay runs the risk of being too specific in that it may fail to detect new or atypical contaminants that might result from some process irregularity or change. For this reason it is recommended that a broadly reactive "generic" or "platform" host cell protein assay be used as part of the final product purity analysis even when a process specific assay is available. If you deem a more process specific assay is necessary, *Cygnus Technologies* is available to apply its proven technologies to develop such antibodies and assays on custom basis.

## Principle of the Procedure

The SF9 Insect Cell Host Cell Protein assay is a two-site immunoenzymetric assay. Samples containing SF9 insect cell proteins are reacted in microtiter strips coated with an affinity purified capture antibody. A second HRP labeled anti-SF9 insect cell antibody is reacted simultaneously, forming a sandwich complex of solid phase antibody-SF9 HCP-enzyme labeled antibody. The microtiter strips are then washed to remove any unbound reactants. After the washes, the substrate tetramethylbenzidine is then reacted. The amount of hydrolyzed substrate is read on a microtiter plate reader and is directly proportional to the concentration of SF9 insect cell HCPs present.

## Reagents & Materials Provided

Component	Product #
<b>Anti-SF9:HRP</b> Affinity purified rabbit antibody conjugated to HRP in a protein matrix with preservative. 1x12mL	<b>F841</b>
<b>Anti-SF9 coated microtiter strips</b> 12x8 well strips in a bag with desiccant	<b>F842*</b>
<b>SF9 HCP Standards</b> SF9 Insect Cell HCPs in a bovine serum albumin matrix with preservative. Standards at 0, 1, 3, 6, 12, 25, 50, and 100ng/mL. 1 mL/vial	<b>F843</b>
<b>Stop Solution</b> 0.5M sulfuric acid. 1x12mL	<b>F006</b>
<b>TMB Substrate</b> 3,3',5,5' Tetramethylbenzidine. 1x12mL	<b>F005</b>
<b>Wash Concentrate (20X)</b> Tris buffered saline with preservative. 1x50mL	<b>F004</b>

\*All components can be purchased separately except # F842.

## Storage & Stability

- All reagents should be stored at 2°C to 8°C for stability until the expiration date printed on the kit.
- Reconstituted wash solution is stable until the expiration date of the kit.
- After prolonged storage, you may notice a salt precipitate and/or yellowing of the wash concentrate. These changes will not impact assay performance. To dissolve the precipitate, mix the wash concentrate thoroughly and dilute as directed in the 'Preparation of Reagents' section.

## Materials & Equipment Required But Not Provided

Microtiter plate reader spectrophotometer with dual wavelength capability at 450 & 650nm. (If your plate reader does not provide dual wavelength analysis you may read at just the 450nm wavelength.)

Pipettors - 50µL and 100µL

Repeating or multichannel pipettor - 100µL

Microtiter plate rotator (400 - 600 rpm)  
Sample Diluent (recommended Cat # 1028)  
Distilled water  
1 liter wash bottle for diluted wash solution

## Precautions

- For Research or Manufacturing use only.
- Stop reagent is 0.5M H<sub>2</sub>SO<sub>4</sub>. Avoid contact with eyes, skin, and clothing. At the concentrations used in this kit, none of the other reagents are believed to be harmful.
- This kit should only be used by qualified technicians.

## Preparation of Reagents

- Bring all reagents to room temperature.
- Dilute wash concentrate to 1 liter in distilled water, label with kit lot and expiration date, and store at 4°C.

## Procedural Notes

1. Complete washing of the plates to remove excess unreacted reagents is essential to good assay reproducibility and sensitivity. We advise against the use of automated or other manually operated vacuum aspiration devices for washing plates as these may result in lower specific absorbances, higher non-specific absorbance, and more variable precision. The manual wash procedure described below generally provides lower backgrounds, higher specific absorbance, and better precision. If duplicate CVs are poor, or if the absorbance of the '0' standard is greater than 0.200, evaluate plate washing procedure for proper performance.

2. High Dose Hook Effect or poor dilutional linearity may be observed in samples with very high concentrations of HCP. High Dose Hook Effect is due to insufficient excess of antibody for very high concentrations of HCPs present in samples upstream in the purification process. Samples with HCP greater than 57 µg/mL may give absorbances less than the 100ng/mL standard. It is also possible for samples to have certain HCPs in concentrations exceeding the amount of antibody for that particular HCP. In such cases, the absorbance of the sample at all dilutions may be lower than the highest standard in the kit, however, these samples will fail to show acceptable dilution linearity/parallelism as evidenced by an apparent increase in dilution corrected HCP concentration with increasing dilution. If a hook effect is possible, samples should also be assayed diluted. If the HCP concentration of the undiluted or less diluted sample is less than a more diluted sample, this may be indicative of the hook effect. Such samples should be diluted at least to the minimum required dilutions (MRDs) as established by your validation

## Assay Protocol

studies using your actual final and in-process drug samples. The MRD is the first dilution at which all subsequent dilutions yield the same HCP value within the statistical limits of assay precision. The HCP value to be reported for such samples is the dilution corrected value at or greater than the established MRD. The diluent used should be compatible with accurate recovery. The preferred diluent is our Cat# I028 available in 100mL, 500mL, or 1 liter bottles. This is the same material used to prepare the kit standards. As the sample is diluted in I028, its matrix begins to approach that of the standards, thus reducing any inaccuracies caused by dilutional artifacts. Other prospective diluents must be tested for non-specific binding and recovery by using them to dilute the 100ng/mL standard, as described in the "Limitations" section below.

## Limitations

- Before relying exclusively on this assay to detect host cell proteins, each laboratory should validate that the kit antibodies and assay procedure yield acceptable specificity, accuracy, and precision. A suggested protocol for this validation can be obtained from our Technical Services Department or our web site.
- The standards used in this assay are comprised of SF9 HCPs solubilized by mechanical disruption and detergent. AAE analysis of the antibodies used in this kit demonstrates that they recognize the majority of distinct PAGE separated bands seen using a sensitive protein staining method like silver stain or colloidal gold. Because the vast majority of HCPs will be conserved among all strains of SF9 this kit should be adequately reactive to HCPs from your strain
- Certain sample matrices may interfere in this assay. The standards used in this kit attempt to simulate typical sample protein and matrices. However, the potential exists that the product itself or other components in the sample matrix may result in either positive or negative interference in this assay. High or low pH, detergents, urea, high salt concentrations, and organic solvents are some of the known interference factors. It is advised to test all sample matrices for interference by diluting the 100ng/mL standard, 1 part to 3 parts of the matrix containing no or very low HCP contaminants. This diluted standard when assayed as an unknown, should give an added HCP value in the range of 20 to 30ng/mL. Consult *Cygnus Technologies* Technical Service Department for advice on how to quantitate the assay in problematic matrices.
- Avoid the assay of samples containing sodium azide (NaN<sub>3</sub>) which will destroy the HRP activity of the conjugate and could result in the underestimation of HCP levels.

- The suggested assay protocol takes approximately 2.5 hours to complete and will yield a sensitivity of <0.5ng/mL. The assay is very robust such that assay variables like incubation times, sample size, and other sequential incubation schemes can be altered to manipulate assay performance for more sensitivity, increased upper analytical range, or reduced sample matrix interference. Before modifying the protocol from what is recommended, you are advised to contact Technical Services for input on the best way to achieve your desired goals.
- The protocol specifies use of an approved microtiter plate shaker or rotator for the immunological steps. These can be purchased from most laboratory supply companies. If you do not have such a device, it is possible to incubate the plate without shaking however, it will be necessary to extend the immunological incubation step in the plate by about one hour in order to achieve comparable results to the shaking protocol. Do not shake during the 30-minute substrate incubation step, as this may result in higher backgrounds and worse precision.
- Bring all reagents to room temperature. Set-up plate spectrophotometer to read dual wavelength at 450nm for the test wavelength and ~650nm for the reference.
- The conjugate will have a cloudy appearance. This is normal and does not indicate contamination. Overtime, you may observe a slight precipitate. This precipitate is inconsequential to assay results. We suggest a simple inversion of the bottle to re-suspend it.
- Thorough washing is essential to proper performance of this assay. Automated plate washing systems or other vacuum aspiration devices are not recommended. The manual method described in the assay protocol is preferred for best precision, sensitivity and accuracy. A more detailed discussion of this procedure can be obtained from our Technical Services Department or on our web site. In addition, a video demonstration of proper plate washing technique is available in the 'Technical Help' section of our web site.
- All standards, controls, and samples should be assayed at least in duplicate.
- Maintain a repetitive timing sequence from well to well for all assay steps to insure that all incubation times are the same for each well.
- Make a work list for each assay to identify the location of each standard, control, and sample.
- It is recommended that your laboratory assay appropriate quality control samples in each run to insure that all reagents and procedures are correct. **You are strongly urged to make controls in your typical sample matrix using**

HCPs derived from your cell line. These controls can be aliquoted into single-use vials and stored frozen for long-term stability.

- Strips should be read within 30 minutes after adding stop solution since color will fade over time.

### Assay Protocol

1. Pipette 100 $\mu$ L of anti-SF9:HRP (#F841) into each well.
2. Pipette 50 $\mu$ L of standards (#F843), controls and samples into wells indicated on work list.
3. Cover & incubate on rotator at 400-600 rpm for 2 hours at room temperature, 24°C  $\pm$  4°C.
4. Dump contents of wells into waste. Blot and gently but firmly tap over absorbent paper to remove most of the residual liquid. Overly aggressive banging of the plate or use of vacuum aspiration devices in an attempt to remove all residual liquid is not necessary and may cause variable dissociation of antibody bound material resulting in lower ODs and worse precision. Fill wells generously to overflowing with diluted wash solution using a squirt bottle or by pipetting in ~350 $\mu$ L. Dump and tap again. Repeat for a total of 4 washes. Wipe off any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step. Do not allow wash solution to remain in wells for longer than a few seconds. Do not allow wells to dry before adding substrate.
5. Pipette 100 $\mu$ L of TMB substrate (#F005).
6. Incubate at room temperature for 30 minutes. DO NOT SHAKE.
7. Pipette 100 $\mu$ L of Stop Solution (#F006).
8. Read absorbance at 450/650nm

### Quality Control

- Precision on duplicate samples should yield average % coefficients of variation of less than 15% for samples in the range of 3-100ng/mL. CVs for samples <3ng/mL may be greater than 15%.
- It is recommended that each laboratory assay appropriate quality control samples in each run to insure that all reagents and procedures are correct.

### Calculation of Results

The standards may be used to construct a standard curve with values reported in ng/mL "total immunoreactive HCP equivalents". This data reduction may be performed through computer methods using curve fitting routines such as point-to-point, cubic spline, or 4 parameter logistic fit. **Do not use linear regression analysis to interpolate values for samples as this may lead to significant inaccuracies!** Data may also be manually reduced by plotting the absorbance values of the standard on the y-axis versus concentration on the x-axis and drawing a smooth point-to-point line.

Absorbances of samples are then interpolated from this standard curve.

### Example Data

Well #	Contents	Abs. at 450nm	Mean Abs.
1A	0ng/mL	0.076	0.078
2A	0ng/mL	0.079	
1B	1ng/mL	0.143	0.141
2B	1ng/mL	0.139	
1C	3ng/mL	0.211	0.208
2C	3ng/mL	0.204	
1D	6ng/mL	0.346	0.337
2D	6ng/mL	0.328	
1E	12ng/mL	0.592	0.587
2E	12ng/mL	0.582	
1F	25ng/mL	1.052	1.030
2F	25ng/mL	1.008	
1G	50ng/mL	1.835	1.780
2G	50ng/mL	1.725	
1H	100ng/mL	3.024	2.982
2H	100ng/mL	2.939	

### Performance Characteristics

Cygnus Technologies has qualified this assay by conventional criteria as indicated below. This qualification is generic in nature and is intended to supplement but not replace certain user and product specific qualification and validation that should be performed by each laboratory. At a minimum each laboratory is urged to perform a spike and recovery study in their sample types. In addition, any of your samples types containing process derived HCPs within or above the analytical range of this assay should be evaluated for dilution linearity to insure that the assay is accurate and has sufficient antibody excess for your particular HCPs. Each laboratory and technician should also demonstrate competency in the assay by performing a precision study similar to that described below. A more detailed discussion of recommended user validation protocols can be obtained by contacting our Technical Services Department or on-line at our web site.

#### Sensitivity

The lower limit of detection (LOD) is defined as that concentration corresponding to a signal two standard deviations above the mean of the zero standard. LOD is 0.3 ng/mL.

The lower limit of quantitation (LOQ) is defined as the lowest concentration, where concentration coefficients of variation (CVs) typically are <20%. The LOQ is 0.5 ng/mL.

## Precision

Both intra (n=20 replicates) and inter-assay (n=10 assays) precision were determined on 3 pools with low (3ng/mL), medium (12ng/mL), and high concentrations (40ng/mL). The % CV is the standard deviation divided by the mean and multiplied by 100.

Pool	Intra assay CV	Inter assay CV
Low	5.5%	4.3%
Medium	3.5%	3.0%
High	4.1%	2.8%

## Recovery/Interference Studies

Various buffer matrices have been evaluated by adding known amounts of SF9 HCPs used to make the standards in this kit. Because this assay is designed to minimize matrix interference most of these buffers yielded acceptable recovery (defined as between 80-120%). In general, extremes in pH (<6.0 and >8.5) as well as certain detergents can cause under-recovery. Organic solvents and high salt concentration can also interfere. In some cases, very high concentrations of the product protein may also cause a negative interference in this assay. Each user should demonstrate that their sample matrices and product itself yield accurate recovery. Such an experiment can be performed by diluting the 100ng/mL standard provided with this kit into the sample in question. For example, we suggest adding 1 part of the 100 ng/mL standard to 3 parts of the test sample. This yields an added spike of 25ng/mL. Any endogenous SF9 HCPs from the sample itself determined prior to spiking and corrected for by the 20% dilution of that sample can be subtracted from the value determined for the spiked sample. The added spike and recovery should be within allowable limits e.g. 80% to 120%. Should you have any problems achieving adequate spike and recovery data, you are strongly urged to contact our Technical Services Department for recommendations on how to overcome sample matrix interference.

## Hook Capacity

Increasing concentrations of HCPs >100 ng/mL were assayed as unknowns. The hook capacity, defined as that concentration that will give an absorbance reading less than the 100 ng/mL standard was >57 µg/mL.

## Ordering Information/ Customer Service

*Cygnus Technologies* also offers kits for the extraction of Host Cell DNA. The following kits are available:

- Residual Host Cell DNA extraction:

Cat # D100W, DNA Extraction Kit in 96 deep well plate  
Cat # D100T, DNA Extraction Kit in microfuge tubes

To place an order or to obtain additional product information contact *Cygnus Technologies*:

[www.cygnustechnologies.com](http://www.cygnustechnologies.com)

Cygnus Technologies, Inc.  
4332 Southport Supply Rd. SE  
Southport, NC 28461 USA  
Tel: 910-454-9442

Fax: 910-454-9443

Email: [techsupport@cygnustechnologies.com](mailto:techsupport@cygnustechnologies.com)





